

iXon EMCCD

Driving the absolute best from
EMCCD technology



iXon

The Industry's Highest Performance Scientific EMCCD Cameras

Andor Technology pioneered the world's first scientific Electron Multiplying CCD (EMCCD) cameras, shipping the initial cameras back in 2000 and winning the Photonics Circle of Excellence award. At that time, Andor coined the name 'Electron Multiplying CCD (EMCCD)', which has been adopted right across this burgeoning industry.

Since then, Andor has consistently set progressively higher EMCCD performance standards with our successive iXon series of deep-cooled, vacuum sealed, quantitative EMCCD cameras. For example, Andor introduced the first back-illuminated EMCCDs in January 2002, alongside our unique 'Baseline Clamp' solution for enhanced quantitative performance. Andor's method for achieving industry-lowest Clock Induced Charge (CIC) was introduced in early 2003 and our benchmark quantitative and linearized EM gain control (RealGain™) and patented EM gain recalibration technology (EMCAL™) was innovated in January 2006.

In 2010, Andor introduced the iXon3 series that brought a number of customer requested features, such as one-click application optimization (OptAcquire) and the ability to calibrate data in either photoelectrons or photons (Count Convert).

The iXon Ultra built upon this rich performance and feature set, raising the bar markedly higher still by driving frame rates up to 3x more than our already industry leading speed performance. This opened new possibilities for the majority of EMCCD-enabled applications that benefit from single photon sensitivity at fast frame rates.


iXon Life, Andor's latest platform addition, is exclusively for fluorescence microscopy applications as is designed to offer EMCCD ultra-sensitivity at a lower price point normally associated with back-illuminated sCMOS cameras, without sacrificing the raw sensitivity and speed required for light starved microscopy applications such as single molecule detection.

“

The fast readout, low noise, and large areas of the iXon cameras allow us to beat the turbulence in the Earth's atmosphere in order to capture the details of the magnetic structuring in the solar chromosphere.

Image courtesy of Mr Kevin Reardon, INAF - Osservatorio Astrofisico di Arcetri, Italy

”



The EMCCD Advantage

Ultimate sensitivity with super fast speeds




Image courtesy of Xiaowei Zhuang, Harvard University, USA

Current trends in photon measurement are placing unprecedented demands on detector technology to perform at significantly higher levels of sensitivity and speed. Electron Multiplying CCD (EMCCD) technology has been designed to respond to this growing need, unlocking new and innovative experimental prospects.

EMCCDs operate by amplification of weak signal events (down to single photons) to a signal level that is well clear of the read noise floor of the camera at any readout speed. Importantly, this ‘on-chip’ amplification process is realized without sacrificing the photon collection capability of the sensor, with back-illuminated sensors offering up to 95% Quantum Efficiency (QE).

iXon Life - The Microscopist’s Choice

In applications such as single molecule microscopy, super-resolution, live cell microscopy (including confocal), calcium signaling, transport/motile imaging and intracellular bioluminescence, weak, rapidly changing fluorescent signals from cells must be dynamically imaged. Andor’s iXon Life platform offers an ideal detection solution. Ultra-sensitive detection capability in fluorescence microscopy facilitates use of lower excitation powers (thereby reducing photobleaching and phototoxicity) and lower dye concentrations.

Since its pioneering introduction in 2000, Andor’s EMCCD technology has been widely and highly successfully employed by microscopists throughout the world, resulting in an outstanding level of representation in high-profile publications.

iXon Ultra - The Physicist’s Choice

The unique high-performance specifications of the optimized iXon Ultra platform have been serving the physical scientist and astronomer in scenarios that demand more than simply an EM sensor in a camera. Andor has worked with numerous scientists to deliver solutions that work for their particular application requirements, such as providing effective charge purging immediately prior to acquisition, specific coatings, coupling to fiber optic scintillators and also specific interface requirements.

As such the Andor iXon brand has been prevalent across a variety of demanding applications, such as photon counting, lucky astronomy, adaptive optics, Bose Einstein condensation (BEC) / ion trapping, single molecule detection / nanotechnology, neutron tomography, X-Ray/Gamma tomography, plasma diagnostics, Raman detection and thermo-luminescence detection.

The iXon Platforms

iXon Ultra

Driving the absolute best from EMCCD technology

The market leading EMCCD camera, Andor’s iXon Ultra has been long recognized as the highest performance, most versatile solution for the most demanding of light starved applications. Available for **both physical and life science** applications, iXon Ultra delivers unbeatable flexibility and functionality, offering **market leading TE cooling, ‘2-in-1’ EMCCD and CCD readout modes** and a host off other high end functionality.

iXon Life **NEW**

Unique price/performance EMCCD exclusively for fluorescence microscopy

Andor’s new **iXon Life EMCCD** platform is available **exclusively for fluorescence microscopy** applications and is engineered to deliver **single photon sensitivity** with absolutely **unparalleled price/performance**. Perfect for single molecule detection and live cell microscopy with minimized phototoxicity or photobleaching, but at a price more normally associated with back-illuminated sCMOS cameras.

-100

Andor’s **UltraVac™** vacuum technology, is critical to ensure both **-100°C deep cooling** and **complete protection of the sensor**.

Model	iXon Ultra 888	iXon Ultra 897
Core Attributes	Field of View, sensitivity and speed	Sensitivity and speed
Sensor Format	1024 x 1024	512 x 512
QE Options	BV, EX2, UVB (Ultra)	BV, EX2, UVB (Ultra)
Pixel Size	13 μm	16 μm
Frame Rate	26 fps (670 fps with 128 x 128 Crop Mode)	56 fps (595 fps with 128 x 128 Crop Mode)
Read Noise	< 1 e ⁻ with EM Gain	< 1 e ⁻ with EM Gain
Minimum TE Cooling Temperature	-95°C	-100°C
CCD Mode (conventional amplifier)	Yes	Yes
Mechanical Shutter	Yes	Yes

Model	iXon Life 888 NEW	iXon Life 897 NEW
Core Attributes	Field of View, sensitivity and speed	Sensitivity and speed
Sensor Format	1024 x 1024	512 x 512
QE Options	BV	BV
Pixel Size	13 μm	16 μm
Frame Rate	26 fps (670 fps with 128 x 128 Crop Mode)	56 fps (595 fps with 128 x 128 Crop Mode)
Read Noise	< 1 e ⁻ with EM Gain	< 1 e ⁻ with EM Gain
Minimum TE Cooling Temperature	-75°C	-80°C
CCD Mode (conventional amplifier)	No	No
Mechanical Shutter	No	No

iXon 888 Models

Field of View and Sensitivity... now 3x faster!



Image courtesy of Jan Liphardt, Stanford University, USA

Available on iXon Ultra and iXon Life platforms, the highly innovative iXon 888 1 Megapixel, back-illuminated EMCCD cameras, offer single photon sensitivity across a large field of view, at 26 fps.

Building on a rich history of first to market innovation, the ‘supercharged’ iXon 888, represents a massive performance boost for the largest available EMCCD sensor.

The iXon 888 has been fundamentally re-engineered to facilitate 3x overclocking of the pixel readout speed to an unprecedented 30 MHz, whilst maintaining quantitative stability, thrusting the full frame rate performance to video rate. Furthermore, Andor’s unique ‘Crop Mode’ can be employed to further boost frame rates from a user defined sub-region, for example pushing a 512 x 512 sub-array to 93 fps and a 128 x 128 area to 670 fps.

With a 1024 x 1024 sensor format and 13 µm pixel size, the resolving power, field of view and unparalleled speed of the iXon 888 render it the most attractive and versatile EMCCD option for demanding applications. These include single molecule detection, super-resolution microscopy, live cell imaging and high time resolution astronomy.



Specification sheets
available at
[andor.com/ixon](https://www.andor.com/ixon)

Key Specifications

Active pixels	1024 x 1024
Pixel size (w x h; µm)	13 x 13
Image area (mm)	13.3 x 13.3
Active area pixel well depth (e)	80,000
Max readout rate (MHz)	30
Frame rates (fps)	26 (full frame) - 9,610
Read noise (e)	< 1 with EM gain
QE max	> 95% (EX2 available)

Features and Benefits

13.3 x 13.3 mm sensor	Largest field of view EMCCD available
Single Photon Sensitive and > 95% QE	Optimal SNR in light starved applications such as single molecule detection and quantum physics
Overclocked to 30 MHz readout	Follow dynamic changing processes.
Crop Mode	Continuous imaging with fastest possible frame rate from centrally positioned ROIs. Highly enabling for live cell super-resolution and much more (e.g. 251 fps with 256 x 256 ROI).
TE cooling to -95°C (deepest cooling with Ultra)	Elimination of dark current detection limit.
Superior Baseline Clamp and EM Stability	Essential for quantitative accuracy of dynamic measurements.
RealGain™	Absolute EMCCD gain selectable directly from a linear and quantitative scale.
Lower Noise CCD Mode (Ultra only)	‘2 in 1’ flexibility. EMCCD for ultra-sensitivity at speed, conventional CCD for longer acquisitions.
EX2 Technology (Ultra option only)	Extended QE response, beyond standard back-illuminated.
Fringe Suppression (Ultra option only)	Reduced etaloning in NIR.
OptAcquire	Optimize the highly flexible iXon for different application requirements at the click of a button.
Count Convert	Quantitatively capture and view data in electrons or incident photons. Applied either in real time or postprocessing, Count Convert does this important conversion for you.
EMCAL™	Patented user-initiated self-recalibration of EM gain.
iCam	Exposure time fast switching provides market leading acquisition efficiency.
Minimal Clock-Induced Charge	Unique pixel clocking parameters, yielding minimized spurious noise floor.
UltraVac™	Critical for sustained vacuum integrity and to maintain unequalled cooling and QE performance, year after year. Seven year vacuum warranty.
Spurious Noise Filter	Intelligent algorithms to filter clock induced charge events from the background. Real time or post-processing.
Direct Data Access (Ultra only)	Camera Link output port to facilitate direct access to data for ‘on the fly’ processing.
Enhanced photon counting modes (Ultra only)	Intuitive single photon counting modes to overcome multiplicative noise. Real time or post-processing.
FPGA Timestamp	Hardware generated timestamp with 10 ns accuracy.


Key Applications

Single Molecule Detection
Cell Motility
Super Resolution (PALM, STORM)
TIRF Microscopy
FRET / FRAP
Fluorescence Correlation Microscopy
Vesicle Trafficking
Microspectroscopy
Lucky Astronomy (Ultra only)
Ion Signaling (Calcium Flux)
Adaptive Optics (Ultra only)
Quantum physics (Ultra only)

The iXon 888 models maintain all the advanced performance attributes and a rich customer requested feature set that have defined the iXon range to date, such as deep vacuum cooling to -95°C, extremely low spurious noise and EM Gain calibration.

Count Convert functionality allows real time data acquisition in units of electrons or incident photons and OptAcquire facilitates one-click optimization of this versatile camera to a variety of application conditions.

Additional features of the iXon Ultra 888 model include a lower noise CCD mode and an additional Camera Link output. The latter offers a unique ability to directly intercept data for ‘on the fly’ processing, ideally suited to applications such as adaptive optics.



Field of View Comparison between iXon 888 and 897 models. The 888 model has a x2.6 greater sensitive area than the 897 model.

iXon 897 Models

Ultimate sensitivity ... supercharged!



Image courtesy of Salvatore Sauro, King's College London, UK

Available on iXon Ultra and iXon Life platforms, the iXon 897 models take the popular back-illuminated 512 x 512 frame transfer sensor and overclock readout to 17 MHz, pushing speed performance to an outstanding 56 fps (full frame), whilst maintaining quantitative stability throughout. Ultimate sensitivity is also attained through deep thermoelectric cooling down to -95°C and industry-lowest clock induced charge noise.

The significant speed boost offered in the iXon Ultra 897 facilitates a new level of temporal resolution to be attained. This is ideal for speed challenged low-light applications such as super-resolution microscopy, single molecule tracking, ion signaling, cell motility, single photon counting, lucky astronomy and adaptive optics.

The extremely low noise of the iXon 897 coupled with the new overclocked speed performance will place this model at the forefront of consideration when it comes to upgrading the high end imaging performance of your laboratory.



Specification sheets
available at
[andor.com/ixon](https://www.andor.com/ixon)

Features and Benefits

Single Photon Sensitive and > 95% QE	Optimal SNR in light starved applications such as single molecule detection and quantum physics
Overclocked to 17 MHz readout	Follow dynamic changing processes.
Crop Mode	Continuous imaging with fastest possible frame rate from centrally positioned ROIs. Highly enabling for live cell super-resolution and much more (e.g. 251 fps with 256 x 256 ROI).
TE cooling to -100°C (deepest cooling with Ultra)	Elimination of dark current detection limit.
Superior Baseline Clamp and EM Stability	Essential for quantitative accuracy of dynamic measurements.
RealGain™	Absolute EMCCD gain selectable directly from a linear and quantitative scale.
Lower Noise CCD Mode (Ultra only)	'2 in 1' flexibility. EMCCD for ultra-sensitivity at speed, conventional CCD for longer acquisitions.
EX2 Technology (Ultra option only)	Extended QE response, beyond standard back-illuminated.
Fringe Suppression (Ultra option only)	Reduced etaloning in NIR.
OptAcquire	Optimize the highly flexible iXon for different application requirements at the click of a button.
Count Convert	Quantitatively capture and view data in electrons or incident photons. Applied either in real time or postprocessing, Count Convert does this important conversion for you.
EMCAL™	Patented user-initiated self-recalibration of EM gain.
iCam	Exposure time fast switching provides market leading acquisition efficiency.
Minimal Clock-Induced Charge	Unique pixel clocking parameters, yielding minimized spurious noise floor.
UltraVac™	Critical for sustained vacuum integrity and to maintain unequalled cooling and QE performance, year after year. Seven year vacuum warranty.
Spurious Noise Filter	Intelligent algorithms to filter clock induced charge events from the background. Real time or post-processing.
Direct Data Access (Ultra only)	Camera Link output port to facilitate direct access to data for 'on the fly' processing.
Enhanced photon counting modes (Ultra only)	Intuitive single photon counting modes to overcome multiplicative noise. Real time or post-processing.
FPGA Timestamp	Hardware generated timestamp with 10 ns accuracy.

Key Specifications

Active pixels	512 x 512
Pixel size (w x h; µm)	16 x 16
Image area (mm)	8.2 x 8.2
Active area pixel well depth (e)	180,000
Max readout rate (MHz)	17
Frame rates (fps)	56 - 11,074
Read noise (e)	< 1 with EM gain
QE max	> 95% (EX2 available)

Key Applications

Single Molecule Detection	Vesicle Trafficking
Cell Motility	Microspectroscopy / Hyperspectral imaging
Super Resolution (PALM, STORM)	Spinning Disk Confocal Microscopy
TIRF Microscopy	Lucky Astronomy (Ultra only)
FRET / FRAP	Ion Signaling (Calcium Flux)
Single Plane Illumination Microscopy (SPIM)	Adaptive Optics (Ultra only)
Fluorescence Correlation Microscopy	Quantum physics (Ultra only)

Performance & Innovations

Industry Fastest Frame Rate

Maximum frame rate performance in EMCCDs is a function of two parameters; (1) Pixel Readout Speed (horizontal); (2) Vertical Clock speed.

The former dictates how rapidly charge is pushed horizontally through the EM gain register and the remaining readout electronics, while the latter dictates the speed at which charge is vertically shifted down through both the exposed sensor area and masked frame transfer area of the chip.

iXon offers industry fastest vertical shift speeds, resulting in faster frame rates and

reduced smearing, significantly faster under commonly employed conditions of sub-array/binning.

Notably, the iXon 897 models overclock the pixel readout speed from the standard 10 MHz to 17 MHz, further boosting the frame rate by >60%, yielding 56 fps (full frame).

The iXon 888 models take this a big step further, thrusting clock speed to 30 MHz! This permits video rate frame rate from this large field of view sensor, and enables as fast as 93 fps from a 512 x 512 ROI in Crop Mode.

Key Features

iXon overclocked up to 30 MHz pixel readout speed: 3x faster full frame rate

Fastest vertical shift speeds yield further speed gains with ROI / binning

Minimized smearing through faster vertical shifts

'Optically Centered Crop Mode' for industry fastest ROI speeds, ideal for live cell super-resolution microscopy

Pushing Frame Rates with Crop Mode




Through this acquisition mode, significant increases in frame rates are accomplished by “fooling” the sensor into thinking it is smaller than it actually is. In standard sub-array/ ROI readout mode each frame still carries the time overhead to readout all pixels to the left and right of the selected area and to vertically shift all pixels above and below the selected area. The charge from these pixels is then dumped before an image is sent from camera to PC. In Crop Mode, the number of pixel readout steps outside of that required to read out the requested sub-array is significantly reduced, resulting in markedly higher frame rates.

However, this mode requires that light is not allowed to fall onto the area of the sensor outside of the defined active sub-area. In optical microscopy, this can be realized in conjunction with the new OptoMask accessory, which inserts easily between the microscope output and the camera. Using the OptoMask, a sub-array can be readily defined through positioning of the masking blades, and a cropped area matched to this in software.

Crop Mode

The active imaging area of the sensor is defined in a way that only a small section of the entire chip is used for imaging.

The remaining area has to be optically masked to prevent light leakage and charge spill-over that would compromise the signal from the imaging area. By cropping the sensor, one achieves faster frame rates because the temporal resolution will be dictated only by the time it requires to read out the small section of the sensor.



Optically Centered Crop Mode: Enabling Live Cell Super-Resolution

The iXon range now comes with ‘Optically Centered Crop Mode’, which gives the user the option to break away from the corner tethered requirement of standard crop mode and select a number of pre-defined ROIs that are located in the centre of the image field.

This is achieved with only minimal sacrifice in achievable frame rate, for example a 128 x 128 optically centered ROI delivering 697 fps. Optically centring of the ROI makes this mode extremely appealing to a number of microscopy techniques, including ‘pointillism’ live cell super-resolution microscopy.

For example, the camera can be operated in full 512 x 512 resolution at a frame rate suited to generation of fixed cell super-resolved images, then Optically Centered Crop Mode can be invoked with a 128 x 128 ROI for generation of super-resolved live cell images showing dynamic events.

Key Features

iXon overclocked up to 30 MHz pixel readout speed: 3x faster full frame rate

Fastest vertical shift speeds yield further speed gains with ROI / binning

Minimized smearing through faster vertical shifts

'Optically Centered Crop Mode' for industry fastest ROI speeds, ideal for live cell super-resolution microscopy

251

The **iXon 888** can achieve a blistering **251 fps** from a **256 x 256 ROI** in **Optically Centered Crop Mode**.

Binning	512 x 512	256 x 256	128 x 128	64 x 64	1024 x 100	1024 x 32	1024 x 1
1 x 1	93 (78)	190 (251)	670 (697)	2,053 (1,319)	259	778	9,690
2 x 2	170 (143)	350 (426)	1,150 (1,019)	3,123 (1,646)	492	1,416	-
4 x 4	291 (245)	601 (653)	1,772 (1,504)	4,109 (1,857)	887	2,370	-

Frame rates achievable by the iXon 888 in Crop Mode - ‘Optically Centered Crop Mode’ frame rates in brackets

Count Convert

iXon offers the capability to quantitatively capture and present data in units of electrons or photons, this important conversion is applied either in real time or as a post-conversion step.

The standard way to present quantitative data in scientific detectors has been in units of ‘counts’, relating to the digitized steps of the Analogue to Digital Converter (ADC) used in the camera. Each Analogue to Digital Unit (ADU) relates to a precise number of ‘photo-electrons’ that were generated originally from photons striking and being captured by the detector pixel.

In the iXon, this conversion factor is very accurately recorded within the camera. Knowing this value, alongside the EM gain (RealGain™) and baseline (bias) offset, facilitates back calculation from the signal

in ADU counts per pixel to the signal in electrons per pixel. Furthermore, knowledge of the Quantum Efficiency (QE) and light throughput properties of the camera at each wavelength enables this process to be taken a step further, allowing the signal to be estimated in photons incident at each pixel, provided the spectral spread of the signal is not too broad.

The Count Convert functionality of the iXon provides the flexibility to acquire data in either electrons or incident photons, with negligible slow down in display rate. Furthermore, the option exists to record the original data in counts and perform this important conversion to either electrons or photons as a post-conversion step, while retaining the original data.

Key Features

Quantify data in electrons or incident photons

Convenient estimate of sample signal intensity at the detector

Real time or post-convert

Reference between different samples, users and set-ups

Meaningful signal relating to PALM/STORM localization accuracy

Performance & Innovations


EX2 Technology - Extended QE from Dual AR sensor coating (iXon Ultra)

iXon Ultra models are now available with a new Dual Anti-Reflection coating applied to the back-illuminated sensor, affording a significant enhancement of the Quantum Efficiency performance.

EX2 technology facilitates broadening of the QE range of the back-illuminated sensors through implementation of a new dual AR coating process, developed by sensor

manufacturer e2v. The net effect is to offer significantly improved sensitivity in both the blue and NIR wavelength regions, whilst maintaining ~ 90% QE across the remainder of the visible region.

EX2 is available on iXon Ultra 897 and 888 models. It is not available on iXon Life models.



Back-illuminated EMCCD sensor QE curves, comparing standard 'BV' mid-band AR coating versus new EX2 dual AR coating. Reference data available at +20°C and -100°C sensor cooling temperatures.

Fringe Suppression Sensors (iXon Ultra)

iXon Ultra models are now available with a new Fringe Suppression property in the sensor design, reducing spatial etaloning effects that can arise through monochromatic imaging in the Near Infra-Red (NIR) wavelength range.

Etaloning is particular to back-illuminated sensors and is caused by interference between reflections off the front and back parallel sensor surfaces. For NIR applications, such as imaging of Bose Einstein Condensates, etaloning effects, often observed more notably beyond ~ 750 nm, can sometimes restrict the ability to perform high integrity quantitative imaging.

Fringe Suppression refers to a sensor design modification that significantly reduces the amplitude of etaloning effects. The design has been implemented by sensor manufacturers e2v and based on a tried and trusted process, successfully validated over several years of CCD manufacture. Fringe Suppression is available on iXon Ultra 897 and 888 models. It is not available on iXon Life models.


“ We have successfully used Andor iXon cameras for many years for super-resolution microscopy, the resolution and sensitivity of these cameras is exceptional. ”

Dr. Mike Heilemann, Institute of Physical and Theoretical Chemistry, Goethe-University, Germany

OptAcquire – Flexibility need not be complicated

The control architecture of the iXon range is extremely flexible, meaning the camera can be adapted and optimized for a wide variety of quantitative experimental requirements, ranging from single photon counting through to slower scan, 16-bit Dynamic Range measurements. However, we are starkly aware that optimizing EMCCD technology is far from trivial, with various set-up parameters influencing and trading off between different camera performance characteristics. We have developed OptAcquire, a unique interface allowing users to conveniently choose from a predetermined list of camera set-up configurations.

The user need only choose how they would like their camera to be optimized, e.g. for 'Sensitivity and Speed', 'Dynamic Range and Speed', 'Time Lapse'. Parameters such as EM gain value, vertical shift speed, vertical clock amplitude, pre-amp sensitivity and horizontal readout speed will then be optimized accordingly, 'behind the scenes'. Furthermore, the option exists to create additional user-defined configurations.



Key Features

- Convenient 'one-click' set-up
- Opens the market leading flexibility of the iXon to less advanced users
- Optimize for range of experimental requirements
- Create additional user defined modes


16-bit

The **iXon Ultra's** flexibility ranges from single photon counting through to slower scan, **high Dynamic Range measurements**.

Performance & Innovations

Deep Thermoelectric Cooling


Single thermal electrons are amplified by the EMCCD gain mechanism. Deep vacuum TE cooling is critical to optimize the sensitivity performance of EMCCD sensors, otherwise the raw sensitivity will be compromised, even under conditions of short exposures.



Images of extremely weak LED signal (signal intensity typical of weak luminescence experiments) acquired with iXon 888 at cooling temperatures -70°C and -95°C (water cooling to achieve latter), 120 second exposure times, sub-region show. The need to push to such deeper cooling temperatures can be readily observed under such extreme low light conditions.

Key Features


- Cooling down to -100°C
- Lowest EM-amplified dark current
- Fewer pixel blemishes (hot pixels)
- Low power consumption vacuum cooling



Line intensity profiles across a row of 512 pixels of the 897 model, taken from DARK images at three different cooling temperatures, 29 ms exposure times. The cleanest noise floor is clearly seen under conditions of deep cooling, even for such short exposure times.

Minimized Clock-Induced Charge

After having minimized dark current through deep cooling, the remaining detection limit in back-illuminated EMCCDs is given by the number of Clock-Induced Charge noise events. Andor's industry-exclusive combination of high resolution clocking parameters and sub-microsecond clock speeds are fundamental to minimizing CIC, enabling truly 'high-end' EMCCD sensitivity to be claimed.




DARK IMAGES taken with the iXon 897 at x1000 gain at different vertical shift speeds, 29 ms exposure time. Cooling temperature was -85°C to ensure minimal dark current contribution.

Spurious Noise Filter

It can still be desirable to optionally filter the remaining spurious noise (Clock-Induced Charge or photons) to give as 'black' a background as possible, eradicating any remaining 'salt and pepper' noise. It is important to utilize noise selection and filter algorithms that are intelligent enough to accomplish this task without impacting the integrity of the signal itself.

This is realized through the new Spurious Noise Filter (SNF) functionality of iXon, which offers the user a choice of advanced algorithms to try. SNF can be applied either in real time or as a post-processing step.



Before and after application of the iXon Spurious Noise Filter


Enhanced Photon Counting (iXon Ultra)

To successfully photon count with EMCCDs, there has to be a significantly higher probability of seeing a 'photon spike' than seeing a dark current/CIC 'noise spike'.

The iXon Ultra 897 combines deepest thermoelectric cooling and low CIC performance, yielding market leading photon counting performance and higher contrast images.

Real-time and post process photon counting...

The advanced photon counting modes of the iXon allow for both real time and post-process photon counting. The latter offers the flexibility to 'trial and error' photon count a pre-recorded kinetic series, trading-off temporal resolution vs SNR.



Performance & Innovations


UltraVac™ Permanent Vacuum Head

It is important that a back-illuminated sensor is housed in a hermetically sealed permanent vacuum head with minimized outgassing, otherwise both cooling performance and the sensor QE will steadily degrade. It is this compelling reason that drove Andor to develop UltraVac.

Andor's proprietary UltraVac process has a proven track record of field reliability, accumulated over more than 15 years of shipping high-end vacuum cameras. UltraVac also enables use of only one input window, improving photon-throughput by 8%.

Key Features

- No QE degradation
- Sustained deep TE cooling
- No maintenance / re-pumping
- One input window
- No condensation



Seven Year Vacuum Warranty

Unlike other vacuum EMCCDs on the market, the iXon family has now been shipping with a vacuum enclosed sensor for almost 10 years, with statistical data that substantiates our extremely robust vacuum claims. With the iXon, Andor are proud to offer an extended 7 year warranty on the vacuum enclosure as standard.

RealGain™, Anti-Ageing and EMCAL™

In early 2006, Andor once again raised the bar by introducing some significant new technology innovations. These particular pioneering steps, were to set new high standards in quantitative EMCCD usage and general EMCCD longevity expectations, which others in the industry are now adopting.

RealGain™

Select absolute EM gain direct from a linear and directly quantitative software scale, x1 to x1000. The EM gain you ask for is the EM gain you get.

Anti-Ageing

Internally configured to significantly inhibit saturation-induced decay of EM gain.

EMCAL™

Innovative user-initiated self-recalibration of EM gain, utilizing a patented method of automated EM gain assessment and Andor's unique Linear and Real Quantitative gain implementation.


Temperature Compensated

Calibration holds across all cooling temperatures. No need to recalibrate on each use in multi-user laboratories and facilities.


Superior Quantitative Stability

The iXon is well regulated in terms of both Baseline (bias offset) rigidity and superior EM gain stability, lending for enhanced quantitative reliability throughout and between measurements.

Baseline Clamp OFF




Baseline Clamp ON



iXon Baseline Clamp (bias stability) in operation

EM Gain Stability



EM Gain stability in the iXon Ultra 897 @ 55 fps. 500 frame kinetic series; frame transfer (overlapped) acquisition; 17.8 ms exposure time; x300 EM Gain.

Direct Data Access (iXon Ultra)

Under standard operation the iXon range uses the USB interface (USB 2.0 or USB 3.0, model dependent) for all control and data transfer with the PC. However, some users require a more direct access to the image data stream, in order that they can perform real-time analysis, possibly using external hardware. Such operation can be particularly important for rapid closed feedback applications such as adaptive optics. Direct real time access to data can also be useful for data intensive applications such as super-resolution microscopy or whole genome sequencing, whereby it can be desirable to carry out real time processing of data on an external GPU, for example.

In order to facilitate such functionality, the iXon Ultra includes an additional Camera Link output port. The Camera Link channel intercepts the image data stream in the camera head immediately after the on-head FPGA processing step, but before the USB frame buffer, therefore undergoes the same amount of on-head image processing.

Key Features

- Direct Data Access via Camera Link output
- Minimal latency or jitter
- USB data stream concurrently accessible
- Compatible with any Camera Link card interface
- On-the-fly data processing
- Ideal for closed loop feedback systems, such as adaptive optics



Performance & Innovations


'2 in 1' Performance - EMCCD and CCD (iXon Ultra)

iXon Ultra 888 and 897 models offer '2 in 1' performance flexibility, in terms of operating as a single photon EMCCD or a low noise conventional CCD, readily user selectable through software selection. Such versatility is attractive in laboratories that can require the camera to operate in low light conditions under both fast and slow frame rates.

In photon starved applications, choosing the EMCCD amplifier usually yields better signal to noise ratio when under faster frame rates conditions (> 1 fps), whereas often the CCD amplifier can yield better signal to noise ratio

when longer exposures can be applied and when the sensor can be read out slowly (i.e. 'seconds per frame' rather than 'frames per second').

This a rule of thumb guide however, and often the choice of amplifier depends ultimately on the light levels available during desired exposure time. Usually it is worth experimenting with the CCD amplifier if the temporal demands are sufficient to readout the sensor at 1 MHz or slower.



Signal to Noise plots, comparing EMCCD mode at any readout speed vs CCD mode at 1 MHz and 80 kHz readout speeds. Higher signal to noise can be secured with the 80 kHz CCD amplifier for light levels greater than 12 photons/pixel, but note that this corresponds to a max frame rate of 0.3 fps.

EMCCD	CCD
✓ Single photon sensitive	✗ 3 to 6 e ⁻ read noise
✗ Multiplication noise	✓ No multiplication noise
✓ Faster frame rates possible	✗ Restricted to slower frame rates

The basic trade-offs between EMCCD and conventional CCD amplifiers

Software Solutions

Andor Solis
Solis is a ready to run Windows package with rich functionality for data acquisition and image analysis/processing. Available on 32-bit and 64-bit versions of Windows (7, 8 and 10). Andor Basic provides macro language control of data acquisition, processing, display and export.

Andor iQ
A comprehensive multi-dimensional imaging software package. Offers tight synchronization of EMCCD with a comprehensive range of microscopy hardware, along with comprehensive rendering and analysis functionality. Modular architecture for best price/performance package on the market.



Andor SDK
A software development kit that allows you to control the Andor range of cameras from your own application. Available as 32 and 64-bit libraries for Windows (7, 8 and 10) and Linux. Compatible with C/C++, C#, Delphi, VB6, VB.NET, LabView and Matlab.

Bitplane Imaris®
Imaris delivers all the necessary functionality for visualization, segmentation and interpretation of multidimensional datasets. By combining speed, precision and intuitive ease-of-use, Imaris provides a complete set of features for handling multi-channel image sets of any size up to 50 gigabytes.

Third party software compatibility
Drivers are available so that the iXon range can be operated through a large variety of third party imaging packages, including:

Metamorph (Molecular Devices Corporation)
NIS Elements (Nikon)
LAS (Leica)
Xcellence (Olympus)
Image Pro (Media Cybernetics)
MicroManager (UCSF)
Till Photonics Live Acquisition (Till Photonics)
Imaging Workbench (Indec)
WinFluor (University of Strathclyde)
Maxim DL (Diffraction Limited)
LabView (National Instruments)
Matlab (MathWorks)



Extensive Imaging Portfolio

The Andor Imaging Range

Have you found what you are looking for? As an alternative to the iXon series, Andor offers an extensive portfolio of high performance low light imaging camera technologies.



iKon CCD Deep cooled, low noise CCD

-100°C cooling
Back-illuminated > 90% QE
1 Megapixel to 4 Megapixel
Enhanced NIR versions
'PV Inspector' model (Optimized for EL / PL in-line inspection)
USB 2.0 true plug and play



iXon EMCCD High performance EMCCD platform

Single photon sensitive and back-illuminated
Industry fastest frame rates
-100°C cooling
Flexible yet intuitive
Quantify in electrons or photons



Zyla sCMOS Fast, sensitive, compact, light sCMOS

< 1 electron read noise @ 30 fps
5.5 and 4.2 Megapixel sensors / 6.5 µm
QE up to 82%
0°C cooling at +27°C ambient
100 fps sustained (10-tap Camera Link)
Cost effective USB 3.0 option
16-bit data range



Neo sCMOS Vacuum cooled, lowest noise sCMOS

1 electron read noise @ 30 fps
5.5 Megapixel / 6.5 µm
-40°C vacuum cooling
30 fps sustained; 100 fps burst
4 GB on head memory
16-bit data range
Fan off vibration free mode

Application & Technical Notes


Andor's high-performance EMCCD cameras have been the technology favorite of the vast majority of EMCCD-enabled laboratories across the globe.

The following section is dedicated to providing a greater depth of understanding of the performance innovations underlying the iXon family of high-end EMCCD cameras, outlining the core technical reasons why Andor is still very much considered the EMCCD industry leaders, notably so in the key areas of sensitivity, speed, stability, longevity, quality and accessibility.

Visit the Learning Center now to discover more at andor.com/learning.

Key iXon Resources include:

- Maximizing frame rate performance in EMCCDs
- UltraVac™ permanent vacuum head and performance longevity
- Count Convert - Quantifying data in Electrons and Photons
- New EMCCD Sensor Enhancements: Extended QE and Fringe Suppression
- Deep Vacuum TE Cooling and Darkcurrent Elimination
- Minimizing Clock Induced Charge - finesse charge clocking
- iXon Ultra Camera Link Output for Direct Data Access
- RealGain™, Anti-Ageing and EMCAL™
- Making Sense of Sensitivity
- iXon Trigger Modes



Notes

[illegible]This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

Customer Support

Andor products are regularly used in critical applications and we can provide a variety of customer support services to maximize the return on your investment and ensure that your product continues to operate at its optimum performance.

Andor has customer support teams located across North America, Asia and Europe, allowing us to provide local technical assistance and advice. Requests for support can be made at any time by contacting our technical support team at andor.com/support.

Andor offers a variety of support under the following format:

- On-site product specialists can assist you with the installation and commissioning of your chosen product
- Training services can be provided on-site or remotely via the Internet
- A testing service to confirm the integrity and optimize the performance of existing equipment in the field is also available on request.

A range of extended warranty packages are available for Andor products giving you the flexibility to choose one appropriate for your needs. These warranties allow you to obtain additional levels of service and include both on-site and remote support options, and may be purchased on a multi-year basis allowing users to fix their support costs over the operating life cycle of the products.



Head Office

7 Millennium Way
Springvale Business Park
Belfast BT12 7AL
Northern Ireland
Tel: +44 (0)28 9023 7126
Fax: +44 (0)28 9031 0792

North America

300 Baker Avenue
Suite 150
Concord, MA 01742
USA
Tel: +1 860-290-9211
Fax: +1 860-290-9566

Japan

5F IS Building
3-32-42 Higashi-Shinagawa
Tokyo 140-0002
Japan
Tel: +81-(0)3-6732-8968
Fax: +81-(0)3-6732-8939

China

Unit 1, Building A,
No. 66 Zhufang Road,
Haidian District,
Beijing 100085
P. R. China
Tel: +86 (0)10-8271-9066
Fax: +86 (0)10-8271-9055

Find us on

