

ZSX Primus IV/IVi

X-RAY FLUORESCENCE SPECTROMETERS

Tube-above/tube-below high-performance models



Rigaku

Leading With Innovation



Industry know-how: **ZSX**

SUPPORTING ANALYSIS USING ZSX GUIDANCE

- Facilitates comprehensive analysis with automatic settings and smart selection features for quantitative applications with minimal operator intervention.

DESIGNED FOR SAFETY


- Tube-above geometry minimizes effects of catastrophic sample failures.


HIGH SPEED AND PRECISION

- High-speed goniometer, high-speed digital counting system.

UNIQUE CAPABILITIES

- Point/Mapping analysis, SQX scatter FP method.






Primus IV features

TUBE-ABOVE EXCITATION SAFELY HANDLES POWDER SAMPLES

- Tube-above optics means the instrument won't be damaged by falling sample particles. Enhanced sensitivity and improved accuracy in the analysis of powder samples made possible by measuring pressed powder samples without using sample films.
- Also enables analysis of liquid samples by using sample cells and liquid sample holders designed for tube above analyses.

Sample spills are prevented with the tube-above optics model

With a tube-above system there is no impact on the optical components due to loose particles falling from a poorly prepared sample. Sample films are often not necessary which can further improve light element sensitivity.



Dust contamination of the vacuum pump is suppressed by the use of a dual vacuum system that equilibrates the vacuum pressure between the sample inlet and the measurement chamber, and by internal filters for powder samples (standard equipment).

SPECIFICATIONS

		ZSX Primus IV	ZSX Primus IVi	
Analysis range ¹		${}_{11}^{\text{Be}} - {}_{96}^{\text{Cm}}$		
Spectral method		Wavelength dispersive		
X-ray generator	X-ray tube	End window type Rh target 4 kW		
	X-ray generator	High-frequency inverter system		
	Heat exchanger	Pure water circulation supplier (built-in)		
Spectrometer	Irradiation method	Tube-above	Tube-below	
	Automatic sample changer	Maximum 48-sample exchange Select 12, 24, 36, or 48 samples (Optional) 96-sample exchange	Maximum 60-sample exchange Select 12, 24, 36, 48, or 60 samples	
	sample loading	Air-lock system		
	Sample size (maximum)	$\phi 52 \text{ mm} \times 30 \text{ mm (H)}$	$\phi 52 \text{ mm} \times 40 \text{ mm (H)}$	
	Primary X-ray beam filter	4 types (Ni400, Ni40, Al125, Al25) (Optional) Be30, for X-ray tube protection		
	Diaphragm ²	(Standard) $\phi 35, 30, 20, 10, 1, 0.5 \text{ mm}$		
	Slit	3 slit exchanger (Standard) Standard and fine (Option) Ultra-light element or ultra-high resolution		
	Goniometer	θ - 2θ independent drive system		
	Continuous scan	0.1° - 600°/min		
	Crystal changer	10 position changer		
	Analyzing crystal	(Standard configuration) LiF (200), GeH, PETH, RX26		
	Optional crystals	LiF (420), LiF (220), RX9, RX4, RX35, RX40, RX45, RX61, RX61F, RX75, RX85		
	Vacuum system	Sample chamber and main chamber-per unit	Sample chamber and main chamber-shared	
		Filter for powder samples		
	Atmosphere	(Option) Automatic helium purge system (Option) Liquid sample holder recognition		
Point/mapping mechanism	r- θ stage (Option) Sample observation system			
Pulse height analyzer	Digital multi-channel analyzer (D-MCA)			
Counting system	Detector	For use with heavy elements	SC (Scintillation counter)	
		For use with light elements	F-PC (Gas flow proportional counter) (Optional) S-PC LE (Gas sealed proportional counter: does not require proportional gas)	

INSTALLATION SPECIFICATIONS

Customer: Depressurization valve

Required power supply	Single or three phase 50/60 Hz Personal computer: 1-phase, 100-240 V, 10 A
Grounding specification	30 Ω or below grounding (Independent)
Cooling water	Temperature: Lower than 30°C Pressure: 0.29 - 0.49 MPa Flow: More than 5 L/min Quality: Equivalent to drinking water
Drained water	Gravity drain
Room temperature	15 - 30°C Daily variation within $\pm 2^\circ\text{C}$
Relative Humidity	10 - 75% RH or less
Gas for detector	P-10 Gas (argon 90% - methane 10% mixed gas) Pressure 0.15 MPa, 7 mL/min * Not required if S-PC LE is selected


Connection port: Tapered female screw for pipes Rc1/4

*1: Depending on crystal configuration
*2: Diaphragm diameters can be selected from the following five combinations

Diaphragm options:
 $\phi 35, 30, 20, 10, 3, 1 \text{ mm}$
 $\phi 35, 30, 20, 10, 3, 0.5 \text{ mm}$
 $\phi 35, 27, 20, 10, 3, 1 \text{ mm}$
 $\phi 35, 27, 20, 10, 3, 0.5 \text{ mm}$
 $\phi 35, 27, 20, 10, 1, 0.5 \text{ mm}$

EXTERNAL SIZE - LAYOUTS Unit: mm

ZSX Primus IV Weight: 620 kg



ZSX Primus IVi Weight: 500 kg

